



MDL



**CrossFire
Commander V6**

Quick Reference Guide

© 1994 – 2002 MDL CrossFire by MDL Information Systems GmbH, Frankfurt am Main, Germany
Beilstein Database: Copyright © 1988-2002, Beilstein-Institut zur Förderung der Chemischen Wissenschaften licensed to Beilstein Chemiedaten und Software GmbH and MDL Information Systems GmbH.
Gmelin Database: Copyright © 2000-2002, Gesellschaft Deutscher Chemiker, licensed to MDL Information Systems GmbH; © 1988-1999, Gmelin Institut fuer Anorganische Chemie und Grenzgebiete der Max-Planck-Gesellschaft zur Förderung der Wissenschaften. All rights reserved.

03/2002

1 GETTING STARTED.....	1		
1.1 LOGON.....	1		
1.2 DATABASE SELECTION.....	5		
1.1.2 SEARCH OPTIONS.....	6		
1.3 EXPERT MODE AND GUIDED MODE.....	7		
2 SEARCH.....	8		
2.1 USING THE TEXT SEARCH.....	8		
2.2 USING THE EASY DATA SEARCH (EDS) OR THE FACT EDITOR	13		
2.3 USING THE STRUCTURE EDITOR.....	19		
2.3.1 REACTION SEARCH BY STRUCTURE.....	21		
2.4 USE A SAVED QUERY.....	23		
2.5 COMBINED SEARCH.....	25		
2.6 START SEARCH.....	26		
3 SEARCH OPTIONS.....	27		
		3.1 QUERY RESULT.....	27
		3.2 AUTOSEARCH.....	28
		3.3 COMBINED SEARCH.....	30
4 DISPLAY HITS.....	31		
4.1 MAIN SCREEN.....	32		
4.2 TREE VIEW.....	35		
4.3 FUNCTION OF THE BUTTONS.....	37		
5. PRINT AND EXPORT.....	40		
5.1 PRINTING.....	40		
5.2 EXPORT AND REPORT.....	41		
5.2.1 <i>One Click Reporting</i>	41		
5.2.2 <i>Reference Reporting</i>	43		
6. DEFAULT SETTINGS.....	44		

1 Getting Started


Welcome to MDL CrossFire Commander V6 SP1. This guide demonstrates briefly the use of the features to create queries and to view the result.

1.1 Logon


To start Commander from whichever platform you are using: double click the icon on your desktop or select in the Windows start menu “**Programs:MDL CrossFire Commander V6**”



The Commander Window looks like this:



- In Commander click the  -Button or select from Menu “**Application:CrossFire**”
- Fill in the login-box as required.



Server User ID and Password

Logon to CrossFire server

Host

User ID

Password

Change password...

Save user id and password


OK Cancel

Click the checkbox if you want your user ID and your password to be saved. You will not be asked anymore for them after clicking the box.

If you need to change your password or enter a new password after you have saved it please proceed the following path:

Select Menu “**Options:Connections..**” and select the profil of which the password is to be changed. Click “**Edit**” and select the checkbox “**Forget user id/password for the server**”. Click “**Ok**”.

You will now be asked for your ID and password upon Logon again.



1.2 Database Selection

	Beilstein Database:	Gmelin Database:
Select Database: <input type="checkbox"/>	<input type="checkbox"/> BS0104AE	<input checked="" type="checkbox"/> GM0103SF

Select Database: These Listboxes select the Databases you will be searching in. Depending on your license you may select one or more (search over multiple database) databases by using the checkboxes. In the case of only one database being checked, the currently active database is displayed bold.

The abbreviation informs you about which version of the database is available:


XX YYYY ZZ

XX:	<i>BS</i> = Beilstein (organic Chemistry)
	<i>GM</i> = Gmelin (inorganic Chemistry)
YYYY:	Year (e.g. <i>01</i>) and update number (e.g. <i>04</i>)
ZZ:	Specification
	Beilstein <i>PR</i> = Plus Reaction
	<i>AB</i> = Abstracts
	<i>AE</i> = EcoPharm
	Gmelin <i>SF</i> = Structure&Facts

1.1.2 Search Options

Text Search: Type in keywords, one or more, to be searched in the basic index.

Example: “malaria”



The screenshot shows a search interface with the following elements:

- Text Search:** Label for the search section.
- in all text fields: Selected search scope.
- in bibliographic fields: Alternative search scope.
- Enter keywords (no field names): Instruction text above the input field.
- YOUR KEYWORD: Placeholder text inside the search input field.
- ...: Ellipsis button to the right of the input field.
- Clear: Button to clear the search input.


Structure Editor: Draw and search for a structure as a molecular species or as component of a reaction.



Structures can be loaded as BCCOM-File or as SD-File using “**File:Open**”, “**File:Execute BCCOM**” or “**File:Execute Query List**”

Enter a **Factual Query** in the Text field, e.g. “bp > 100” or “Author is Maier”

Use the **Fact Editor (Table)**, the **Easy Data Search (EDS)** or the predefined forms to conduct a search for specific facts.



 Use a combination of these interfaces for improved search results

1.3 Expert Mode and Guided Mode

This guide describes the use of Commander in “Guided Mode” in which the Client interacts with the user to ensure a query according to the wish of the user.

If you do not want this interaction you can switch to “Expert Mode” by checking the corresponding checkbox. In this mode whatever was typed in is used directly as the query, it is designed for use by users with highly advanced knowledge of the CrossFire System.




2 Search

2.1 Using the Text Search

The Text Search is designed to search a keyword or expression against the basic index of the database. The result provides a general overview to a topic.

Select the level of your search:



The diagram shows a grey box containing two radio button options: "in all text fields" (selected) and "in bibliographic fields". Two arrows point from these options to callout boxes. The first arrow points from "in all text fields" to a box containing the text "Search your keyword in all text fields". The second arrow points from "in bibliographic fields" to a box containing the text "Search your keyword in the citation index only".

The level determines the indexes which will be browsed with the text you enter in the query field. The Basic Index of the databases is searched for the words entered in the field.

Type your Keywords into the field, you can enter one or more words:




The screenshot shows the "Text Search" interface. On the left, the text "Text Search:" is followed by two radio button options: "in all text fields" (selected) and "in bibliographic fields". To the right of these options is a text input field with the placeholder text "Enter keywords (no field names):". The input field contains the text "Aldol Condensation". To the right of the input field is a small button with three dots "...". To the far right is a "Clear" button.

Examples: "Malaria" or "Antimalaria", "Aldol Condensation", "Cancer Treatment", "Effect of Melatonin on sleep disorder". We recommend to limit the search term to 5-6 words maximum.

Submit the query by carriage return or F7 or press “Start Search”. If more than one word is typed an AutoSearch search will be performed.

If you have selected “in bibliographic fields” only the Basic Index Citations will be searched. The result will be a list of references.

If you have selected “All Text Fields” you will be prompted to determine the context you want to apply for your search.



Your Query → Your query: Aldol Condensation

Result in Database: BS = Beilstein, GM= Gmelin
Check the Radiobutton to see the results in the different databases

Listing word by word
single terms can be selected/deselected

Occurrence of Keywords in Index
Substance / Reactions / Citations

Words:	Substances	Reactions*	Citations
aldol	58	149	3508
condensation	806	13411	14614


Options: Set Truncation and Boolean Operators

Options... AutoSearches: near
Truncation: left

Help Cancel Start Search

Choose “Substance”, “Reactions” or “Citations” to determine the context you want to be searched. The number of occurrences is shown for each word in the column of the context for each word separately. Choose the context you want by clicking the header.

Click the Options button to change the default settings:



Truncation:

Truncation can be set off to search for an exact term.


A star at the end of the keyword can be used, it will switch on the option “Use right Truncation” automatically.

Note: Right Truncation will automatically suppress the word endings “-s” (for plural) and “-ing” to ensure the query regarding the structure of our databases.

Left Truncation will lead to a “> number” in the Information Box, correct numbers can’t be displayed for technical reasons

- For more sophisticated search Boolean operators can be used.

If a “query-like” text is introduced, for example “bp=100”, the user will be prompted:



If “Yes” is chosen Commander translates the Input into the Fact Editor Field. Please press “Start Search” again to submit your Query.

If “No” is chosen Commander will execute the entered text as a text search.

2.2 Using the Easy Data Search (EDS) or the Fact Editor

The upper field of the EDS is for expert search using the Commander's operators and field names:

EDS (Easy Data Search):
 Beilstein Database
 Your query or field name:


bp>150

It is possible to enter the query as a sentence like “Author is Mayer”, “Effect is anticancer”.

EDS (Easy Data Search):
 Beilstein Database
 Your query or field name:

effect is anticancer

In this case Commander translates the sentence into the query language and will prompt the user if more than one field could be applied:



Show form to formulate your query
(see below)

Show all database fields
Use if suggestions do not match your intention

Possible Fields
select one or double click to open EDS or table

Field Values found containing terms of your query. Double Click to start search

Field Search in database Beilstein

Your query:
effect is anticancer

Show Form... Show All Fields...

Open an EDS form

PHARM.E: Effect
PHARM.EP: Endpoint of Effect
ECT.E: Effect
ECT.EP: Endpoint of Effect
OPT.KW: Description
Magnetic birefringence (Cotton-Mouton effect)
Electric birefringence (Kerr effect)
Cotton Effect (abnormal curve)
ESR.KW: Description
1H-electron Overhauser effect
ROT.KW: Description
Stark effect
RAMAN.KW: Description
Raman resonance effect
OSM.KW: Description
M... effect

Help Cancel Search Now

The button **Show Form...** transfers the selected field into a table for further refining. If “Open an EDS Form” is checked an EDS Form is used if available.

Add additional search term

Look up possible values

Bioactivity

EcoPharm

Find all compounds with Pharmacological and Toxicological data

or specify the details you are interested in:

Effect: starts with

Species: starts with

Route of application: starts with

Type (IC50, LD50): starts with

Further Details

Keyword: starts with

Translation into EDS Form

Field Name

Qualifier

Value
(your search term)

Searching for a database field

Field name:

Field content: starts with

(press button for more fields) (a value is required)

Further Fields/Parameters:

Add additional search term


Delete Line

Look up possible values

Translation into Table

A combination of sentences can be used as well, for example “effect is anticancer and species is mouse”. However, due to the complex structures of our databases it is not recommended.

Deselect “Open an EDS Form” and click “Show Form”.




Successful Translation of a two sentence query

Note: Operators to connect the two sentences will be set by Commander according to the nature of the query. They can’t be changed.


Use Predefined Forms:

The lower field enables you to use predefined forms to start a factual search. Double Click or select and click “Open” to enter data:



or double click “Fact Editor (Table) to search for specific values of specific Index Fields:

The Fact Editor Table displays the data structure of the selected Database and enables direct lookup of values for a field:



Operators: and, near, next, not, or, proximity

Double Click (or press F2) to search for a field

Double Click (or press F2) to search for a value

Operator	Field Name	Field Value (press F2 for more)
	pharm.e	anticancer
+or	pharm.e	cancer antigens
and		

Available Fields

- Beilstein Datastructure EcoPharm
 - Basic Indexes (BI)
 - Substance (IDE)
 - Bibliographic Information (BIB)
 - References (REFX)
 - Chemical Data (CHE)
 - Physical Data (PHY)
 - Pharmacological and Ecological Data (PED)
 - Pharmacological Data (PHARM)**
 - Effect (PHARM.E)
 - Endpoint of Effect (PHARM.EP)
 - Species or Test-System (PHARM.SP)
 - Sex (PHARM.S)
 - Route of Application (PHARM.BA)
 - Concentration (PHARM.C)
 - Kind of Dosing (PHARM.KD)
 - Exposure Period (PHARM.FX)

Expand field Effect

anticancer


Occurrence	Content
64	anticalmodulin activity
2	anticancer
2	anticancer, alkylating
15	anticancer, antibiotic
1	anticancer, antibody
21	anticancer, antimetabolite
7	anticancer, formulation
2	anticancer, hormonal
2	anticancer, immunological
386	anticancer, other
2	anticancer, photosensitizer
100	anticancer, radio/chemoprotective
22	anticancer, radio/chemosensitizer
2	anticarcinogenic
25	anticatalensic

2.3 Using the Structure Editor

Two Structure Editors are available in this version: ISIS/Draw 2.4 (default) and the MDL CrossFire Structure Editor. You can select the tool of your choice at Menu “**Options: Structure Editors...**”

Press **Structure Editor:**  or double click in the *Structure Editor Window* to open the Drawing tool, and draw the structure or fragment, or the reaction or parts of a reaction you want to search for.

Go back to Commander by clicking on  (MDL ISIS Draw) or  (MDL CrossFire Structure Editor). Once a structure is loaded in Commander general options for the search can be set on the right side of the screen



Search:

as structure

as reactant

as product

Select the role of the drawn structure:

- *as structure to retrieve data about the compound*
- *as reactant to view the reactivity of the compound*
- *as product to retrieve preparations of the compound*

Stereo:

off

off

absolute

relative

racemic

Set stereochemistry options

Free Sites:

hetero atoms

all atoms

Allow free sites on hetero atoms or all atoms.

No selection: exact search

Allow:

addl. fragments

isotopes

charges

radicals

additional rings


Use checkboxes to select the options you want to apply .

To toggle all options on/off hold CTRL and click one of the options


2.3.1 Reaction Search by Structure

Reactions are typically searched by structure. Full reactions, when starting material and product is known, can be searched by drawing the reaction in your preferred structure editor (see xxx)

Reaction Drawing in ISIS/Draw





Reaction Drawing in CrossFire Structure Editor



Commander recognizes automatically the role of the structure as **Reactant** and **Product** by the position of the arrow.

Half reactions can be drawn in similar way, however the structure options window described above can be used for more flexible query formulation:


Draw a structure and load it into Commander. Select the role of the structure using the radio buttons.

2.4 Use a saved query


A query including all options can be saved for future use or as template using “**File:Save**” or “**File:Save as..**”. from Commander. These can be text, fact or structure queries.

Using “**File:Open**” these can be loaded again into Commander. In case of a structure the options for the saved file can be different from the current settings. If this is the case Commander prompts you to select the options you want to apply:



The behavior of this prompt can be customized:

Use “**Query:File/Open...Options**” to select



2.5 Combined search

MDL CrossFire Commander V6 offers the unique opportunity cross search multiple Databases. Depending on your license you can choose MDL's Beilstein and Gmelin Databases; if available customer Databases can be included as well.

Click the checkbox of each Database to select or deselect it.

	Beilstein Database:	Gmelin Database:
Select Database:	<input checked="" type="checkbox"/> BS0104AE	<input checked="" type="checkbox"/> GM0103SF

If more than one Database is selected the EDS are not active as they are specific for each single database. Therefore a multiple search can only be carried out with a structure and/or a text as query.


The results from each database will be displayed in parallel, in separate windows on a single screen, since different fields are connected with the content of each database.

2.6 Start Search

☞ Searches can only be started from Commander.

- use the “**Start Search**” Button
- Hit “Enter” on your keyboard
- Press F7

to start your search.




3 Search Options

Depending on the search option you used, Commander displays the result as described below. You always find a Button **Change Query** or **Back to Query** that enables you to rework your query when the result is insufficient.


3.1 Query Result

The result of the search without further options is displayed like:



3.2 AutoSearch

If “AutoSearch” is enabled Commander will refine the structure or the text entered as the query step by step to find the closest hit. “AutoSearch” can be enabled/disabled in the Menu “**Query**”.



AutoSearch can be stopped after Commander found the first Hit. If this option is not selected Commander will use all possible variations of the query.




Variations when using a structure:

original query – free sites on hetero atoms – free sites on all atoms – allow additional rings

Variations when using a text search (more than one word):

Boolean Operators “near” – “proximity” – “and”

The result is displayed, including information on the variation of the query

Select	Hits	Note	Query
01 <input type="checkbox"/>	0	Original Query	
02 <input checked="" type="checkbox"/>	5	Relax structural restraints: free sites on hetero atoms	
03 <input checked="" type="checkbox"/>	5	Relax structural restraints: free sites on all atoms	impl. free sites 

Information:
Variation of the Query

Hitset Number
Select Box
Number of Hits

Query used

Use the Checkboxes 02 to select one or more hitsets for display.

3.3 Combined Search


The result of a search in multiple Databases (combined search, see 1.1.1) is displayed in the same fashion.

Select	Hits	Note	Query
01 <input checked="" type="checkbox"/>	415	Beilstein:	Factual Query: "malaria* or protozoa*"
02 <input checked="" type="checkbox"/>	7	Gmelin:	Factual Query: "malaria* or protozoa*"

Use the checkboxes to select the Hitsets for display.

4 Display Hits

From the display of the query result click  or  to open “Display Hits”.







The screenshot shows the 'Display Hits' application window. The title bar reads 'Display Hits - [Beilstein:Q13 hit 1 of 5]'. The menu bar includes File, Edit, Task, View, Options, Window, and Help. The toolbar contains various icons, with the 'View all' button highlighted by a red box. On the left, a 'Tree View' pane shows a hierarchical list of search results, including 'Gmelin(GM0103SF)', 'Beilstein(BS0104PR)', and 'Beilstein(BS0104AE)'. The main display area is divided into four quadrants, each showing a chemical structure and its corresponding hit number and formula: Hit 1: 8687490 C₁₀H₉N₃O₂, Hit 2: 8699968 C₁₄H₁₉N₃O₂, Hit 3: 8705408 C₁₆H₁₇N₃O₃, and Hit 4: 8711363 C₁₇H₁₃N₃O₃. A red 'X' icon is visible in the top toolbar, and a 'View all' button is also present. The status bar at the bottom shows 'idle', 'STRUCT', and 'Subs'.


4.1 Main Screen

The main screen starts on default with the Short Display of the results. Depending on the query the hits are automatically displayed as

Structure

	BRN	Molecular Formula	
Hit 1:8687490 C ₁₀ H ₉ N ₃ O ₂			
			
Hit 2:8699968 C ₁₄ H ₁₉ N ₃ O ₂			Note: If these are displayed in red: Biological Data are available!
			Structure
			
Hit 3:8705408 C ₁₅ H ₁₇ N ₃ O ₃			
			
Hit 4:8711363 C ₁₇ H ₁₃ N ₃ O ₃			
			

as Reaction



-or as citations

Hit 1:5500578 Journal; Pessi, Antonello; Bianchi, Elisabetta; Bonelli, Fabio; Chiappinelli, Lorella; JCCCAT; J.Chem.Soc.Chem.Commun.; EN; 1; 1990; 8-9.

Hit 2:5501252 Journal; Escande, C.; Chevalier, P.; Verdier, F.; Bourdon, R.; JPMSAE; J.Pharm.Sci.; EN; 79; 1; 1990; 23-27.

Hit 3:5503456 Journal; Narita, Mitsuaki; Takegahara, Hitoshi; Ono, Seiji; Sato, Hisaya; BCSJA8; Bull.Chem.Soc.Jpn.; EN; 63; 2; 1990; 484-488.

Hit 4:5509164 Journal; Jung, Mankil; Li, Xun; Bustos, Daniel A.; ElSohly, Hala N.; McChesney, James D.; Milhous, Wilbur K.; JMCMAR; J.Med.Chem.; EN; 33; 5; 1990; 1516-1518.


Hit 5:5525841 Journal; Mansor, S. M.; Ward, S. A.; Edwards, G.; Hoaksey, P. E.; Breckenridge, A. M.; JPPMAB; J.Pharm.Pharmacol.; EN; 42; 1990; 428-432.

Hit 6:5528096 Journal; Hooper, M; Kirby, GC; Kulkarni, MM; Kulkarni, SN; Nagasampagi, BA; et al.; EJMCA5; Eur.J.Med.Chem.Chim.Ther.; EN; 25; 9; 1990; 717-723.



Double Click the item you want to see or select it and click  for details.

4.2 Tree View



On the left side of the screen the tree view provides fast navigation through saved Hitsets and Session Hitsets.

Select the Hitset you want to see and the main screen switches to Short Display.


Click on  to drill down to details, ultimately to field level.

Click on these details to display them in the main screen.

The Tree View can be used to change the context of a Hitset from the original to Citations, Reactions or Substances. Do a right click onto the Hitset and select “**Convert...**”



Choose the extent of the conversion and the context, which you want to put the hitset in (here converting a hitset in Citation Context to Substances or Reactions)



4.3 Function of the Buttons

The buttons in the Display Hits window have the following functions:



Export selected data. Select from predefined Export definition or create your own (choose Settings in Drop down menu)



Goto First / Previous / Next / Last hit of the selected hitset



Go to Commander. You can also press F7



Switch on/off Display of structure in Main Screen



Copy structure to Structure Editor (opens the Editor)



Toggle on/off Short Display



Click to view the search query used




Define the view on the data:

- Identification (Substance/Reaction/Reference and Field Availability List)
- All Fields
- Hit Only (Substance/ Reaction/Reference and Abstract)
- Define your own View File

To define your own view click on Button “View” and select “Define...”. In the window which opens you will find your saved views. Select the View File and press “Apply”. Press “Edit” if you want to Edit the View File. Press “Remove” if you want to delete the View File.

To create a new one press “New...” Select the Field you want to apply in the left screen and press “Copy” to add it to the View File.



5. Print and Export

5.1 Printing



Choose “File:Print...” or [CTRL+P] or press  and select the range you want to Print (Current Hit; Hit 1 – X). Select your options and the layout and press OK. Commander will print the View of the Main Screen (see 3.3.3) on your standard printer.

If you want to print only some facts from a hit, mark the desired detail at the check box on the top right

Reaction 1 of 2			<input checked="" type="checkbox"/>	
Reaction ID	8713044			
Reactant BRN	8687490 1-hydroxy-2,3-dihydro-1 <i>H</i> -3,4a,5-triaza-fluoren-4-one			
Product BRN	8682045 3 <i>H</i> -3,4a,5-triaza-fluoren-4-one			



Select one or more facts and start printing as described above.

5.2 Export and Report

All information in your hits can be exported to other applications such as Microsoft Word, Microsoft Excel or your preferred HTML-Browser. You can export full hitsets, a range of hits or extract information from a hit.

Select the information and use the menu “**File:Export**” or the Button  to define the format for export.

You will be prompted for the settings of your export. Choose one of the predefined settings for Export in various formats or “Settings..” to define your own.

The selected export destination will be opened with the export file when finished.


5.2.1 One Click Reporting

In the Full Display view, each field is equipped with a Title Bar on top that is used to select factual data (see 4.1) and to perform “One-Click-Reporting”:

Button One-Click-Reporting

Reaction 1 of 2		  
Reaction ID	8713044	
Reactant BRN	8687490 1-hydroxy-2,3-dihydro-1 <i>H</i> -3,4a,5-triaza-fluoren-4-one	
Product BRN	8682045 3 <i>H</i> -3,4a,5-triaza-fluoren-4-one	

Clicking on the button or a right mouse click anywhere inside the fact opens a menu:



The selected fact will be added at the end of the defined report (see first row). Use the options to operate the Report function.

To view the report select “Open the Report”. To print the selected facts choose “Print” from the “File”-menu and check the range “Selected Facts”.

5.2.2 Reference Reporting

MDL Commander V6 offers the possibility to export reference lists in various ways:

- If a hit is in Substance context the last item of the tree of the substance is “Reference”. Select it and a deduplicated list of all references for this hit are displayed in the Main display. You can report, export and print it in the way described earlier.
- If you have a hitset in citation context you find a predefined export setting (see Chapter 5.2) for export of references to Endnote, Procite or other reference managing applications.
- From each hitset you can export a list of all references as HTML. You will find a predefined export setting if you click the “Export”-Button. A framed document is produced , in which you find the compounds you exported on the left side, select them by clicking on the hit number. The facts of this hit will be displayed in the right frame. On top of the frame you find a hyperlink “Reference List”. Clicking this link leads to a deduplicated list of all the references of the exported compounds.

6. Default Settings

Commander

General:

Application:	No Application
AutoSearch	Enabled
Stop AutoSearch after 1 st Hit	Disabled
Save Password	Disabled
Startup Task	No Application
Expert Mode	Off
Save Settings on Exit	Enabled
Save History on Exit	Enabled
Automatic Context Selection	Enabled
Show Box before a multiple Database Search	Enabled
Show Startup Screen	Enabled
File/Open...Options	Show query options select box

Structure:

Structure Editor

MDL ISIS Draw

Options:

Search

As structure

Stereo

Off

Free Sites: hetero atoms

Disabled

Free Sites: all atoms

Disabled

Allow addl. Fragments

Enabled

Allow isotopes

Enabled

Allow charges

Enabled

Allow radicals

Enabled

Allow additional rings

Enabled

Define Structure and Reaction View:

Mean Bond Length

12 mm

Atom View

Atom Symbol

Disabled

Without C

Enabled

C as Dot

Disabled

Atom Number

Disabled

Atom Attributes

Number

Disabled



Mass	Enabled
Topology	Enabled
Charge	Enabled
Radical	Enabled
H-Count	Enabled
Free Sites	Enabled
Valency	Enabled
Attachement Points	Enabled
Attachement Number	Enabled
Stereo	Disabled
<i>Bond Attributes</i>	
Topology	Enabled
Stereo	Enabled
Site Number	Enabled
<i>Reaction Attributes</i>	
Atom Mapping	Enabled
Reaction Center	Disabled
New Bond Type	Enabled
 <i>Text Search:</i>	
Text Search	In all text fields

Options:

Words are close to each other	Enabled
Words are in the same fact	Enabled
Words are in the same document	Enabled
Use right truncation	Enabled
Use left truncation	Disabled

Display Hits

Short View	Enabled
Tree View	Enabled
Open Hitset in Short Display View	Enabled All
Save Settings on Exit	Enabled
Include Structures	Disabled
Include Field Availability	Enable
Selected Hits only	Disabled
Status Bar	Enabled

Distributed by:



MDL Information Systems GmbH
Theodor-Heuss-Allee 108
D-60486 Frankfurt
Germany

Headquarters:

MDL Information Systems Inc.
14600 Catalina Street
San Leandro, CA 94577
USA

Phone: +1-800-275 6094 (toll free)
Fax: +1-800-223 0801 (toll free)
e-mail: info@mdl.com

Europe Customer Support:

Phone: +41 61 486-8880
Fax: +41 61 486-8889
e-mail: hotline@mdl.com

US Customer Support:

Phone: +1-800-326 3002 (toll free)
+1-510-895 2213 (outside USA)
Fax: +1-510-895 6092 (outside USA)
e-mail: techsupp@mdl.com

UK Customer Support:

Phone: +44 1276 701-520
Fax: +44 1276 701-501
e-mail: Ukhotline@mdl.com